Analyzing the Polarization State of Light through the Fourier Series

Alamdar Hussain and Muhammad Sabieh Anwar

Syed Babar Ali School of Science and Engineering, LUMS

Version 1; May 20, 2014

One can completely determine the polarization of light simply by using a polarizer and a quarter wave plate (QWP). In this experiment you will develop a method to generate as well as analyze different polarization states of light. Furthermore you will learn to utilize matrix manipulations for solving systems of equations.

Essential pre-lab reading: "*Introduction to Optics*" by F. L. Pedrotti, L. S. Pedrotti and L. M. Pedrotti, Pearson Education, 2008; (Chapter 14).

"Linear Algebra and its Applications" by Gilbert Strang, Cengage Learning, 2006; (Section 3.3, Page 160 to 168).

1 Test Your Understanding

1. Formulate a method to generate circular and elliptical polarizations.




Figure 1: The polarizer and analyzer are oriented along the same direction. The retarders QWP 1 and QWP 2 are at α and β with respect to the first polarizer respectively. The red line shows the perceived path of laser.

- 2. Using Jones calculus, calculate the intensity output from the setup shown in Figure 1.
- 3. Simplify and write the answer in the form of a Fourier series of β . The Fourier series will have the form

$$f(\beta) = C_0 + \sum_{n=1}^{\infty} C_n \cos(n\beta) + \sum_{n=1}^{\infty} S_n \sin(n\beta).$$

What are the Fourier coefficients? Your answer should show that only four of the Fourier coefficients are non-zero.

- 4. Using matrices devise a method for solving a system of equations with more number of equations than the number of unknowns.
- 5. Devise a strategy for finding α using β . Finding α is equivalent to finding the polarization state of light generated by the first set of QWP and the polarizer.

2 The Experiment

A HeNe laser beam is incident on a polarizer. The polarized light then passes through two QWPs, a neutral density filter and an analyzer. Finally it is received by a photodiode.

Orient the polarizer such that it produces a horizontal polarization. All angles will then be measured with respect to this polarizer. Generate a polarization by rotating the first QWP to some angle α between 0° to 45° (try choosing angles away from these limits). Then change the angle β of the second QWP from 0° to 360° in steps of 20° and record the photodiode output for each step. If β is increased in clockwise direction then all choices of α in the clockwise direction are positive and anticlockwise ones are negative.

For finding α , use the expression for output intensity in the form of Fourier series of β and find the coefficients of this series using the pseudo inverse matrix technique (e.g. refer to the appendix and the pre-lab reading). All four coefficients are functions of α only. Use them to find α and hence the polarization state of light.

Q 1. Which coefficients can give information about the handedness (*right* or *left*) of the polarization?

Q 2. Does the calculated value for α match with the experimentally selected one?

Q 3. Repeat the experiment for three different choices of α .

Appendix

Suppose we are given equal number of equations as there are unknowns then the solution is simply given by

$$A_{m \times m} X_{m \times 1} = B_{m \times 1}$$
$$X_{m \times 1} = A_{m \times m}^{-1} B_{m \times 1}$$

where A is the matrix of coefficients, X is a column vector of unknowns or variables and B represents the right hand sides of the equations.

However if there are more number of equations than unknowns then A becomes a rectangular matrix and A^{-1} becomes non trivial. Suppose there are m equations with n unknowns such that m > n then we have

$$A_{m \times n} X_{n \times 1} = B_{m \times 1}$$
$$A_{n \times m}^T A_{m \times n} X_{n \times 1} = A_{n \times m}^T B_{m \times 1}$$
$$X_{n \times 1} = (A_{n \times m}^T A_{m \times n})^{-1} A_{n \times m}^T B_{m \times 1},$$

where $(A^T A)^{-1} A^T$ is called pseudo inverse of A.