

Primer on the Lock-in Amplifier*

Muhammad Sabieh Anwar
LUMS School of Science and Engineering
Version 2016-v1

December 6, 2016

The resistance of a copper wire is so small that it cannot be accurately and reliably measured by ordinary voltmeter. The same is true for a gold nanowire tethered between two supports. In experiments, we repeatedly come across situations where we want to measure the response of a system and the response is feeble or is buried in excessive noise. Through this experiment we introduce the art of precise measurements using a lock-in amplifier which is an integral component of several experiments in the physics lab. As such, the present experiment is a pre-requisite to some of the advanced experiments in this course.

KEYWORDS

Lock-in amplifier · White noise · Flicker noise · Thermal noise · Phase sensitive detection · Low-pass filter · Optical chopper

SCHEDULED TIME 1 WEEK

1 List of Equipment

1. Lock-in Amplifier, SR510, SRS
2. Oscilloscope
3. Signal Generator
4. Electrical Breadboard
5. BNC cables
6. Connecting Wires
7. Resistors
8. Optical Chopper, SR540, Stanford Research Systems (SRS) (Optional)
9. HeNe Laser with Mount, HR020, Thor labs (Optional)

*No part of this document can be reproduced or re-used without the explicit permission of Dr. Muhammad Sabieh Anwar. If you like to use this experiment, please contact him at sabieh@lums.edu.pk.

10. Photodetector, 818-SL, Newport (Optional)
11. Polarizers, LPVIS050, Thor Labs (Optional)
12. Optical Breadboard (Optional)

2 Objectives

In this experiment, we will,

1. eliminate unwanted noise, extract signals buried in excessive noise using lock-in amplifier.
2. learn about different types of noise, and see how noise can obscure the required signal and how phase-sensitive techniques can be used to recover the signal, and
3. find the minute resistance of a conducting wire.
4. optionally, we will use an optical chopper to measure weak optical signals modulated at a certain frequency.

References

- [1] K.G. Libbercht, E.D. Black and C.M. Hirata, "A basic lock-in amplifier experiment for the undergraduate laboratory", Amer. J. Phys. **71**, 1208-1213 (2003).
- [2] Paul Horowitz, W. Hill *The Art of Electronics, 2nd Edition* {Cambridge}, Pg. 184,430-433.
- [3] Manual, Lock-in Amplifier, Model SR510, Stanford Research Systems.

3 Theoretical Introduction

3.1 Noise and its types

Noise obscures the signal of interest [2].

White noise has all frequency components, ranging from zero frequency (DC) to infinite frequencies. The graph between the noise power density versus frequency would be a constant line. Thermal or Johnson noise is a type of white noise. It arises from the thermal fluctuations in a resistor at finite temperature. The rms amplitude of thermal noise voltage is:

$$V_{rms,thermal} = \sqrt{4k_BTRB},$$

where R is the resistance and B is the bandwidth of measurement, T is the absolute temperature, and k_B is the Boltzmann Constant.

Shot noise is another type of white noise which arises due to fluctuation of single charge carriers, e.g. electron flowing through a single electron transistor or a highly attenuated beam of practically single photons impinging on a photodetector.

Pink noise has a power spectral density that decreases with frequency. It has a frequency dependence of $1/f$. It is sometimes also referred to as flicker noise.

“Intrinsic” noise sources include resistors, vacuum diodes, p-n junctions etc. It is hard to eliminate these services. In addition to the intrinsic noise sources, there are a variety of “non-essential” noise sources, i.e., those noise sources which can be minimized with good laboratory practices. Some of these extrinsic noise sources are highlighted below.

Capacitive Coupling A voltage on a nearby piece of apparatus can couple to a detector via stray capacitance. Although the stray capacitance maybe small, but when coupled in, the resultant noise may be larger than weak experimental signal. It can be cured by removing or turning off the interfering noise source, or by installing capacitive shielding which can be done by placing the setup in a metal box. Consider Figure 1(a).




Figure 1: (a) Diagrammatic representation of capacitive noise coupling, (b) inductive noise coupling, and (c) resistive noise coupling.

Inductive Coupling Extrinsic noise can also couple to the experiment via a magnetic field. A changing current in a nearby circuit gives rise to a changing magnetic field which induces an emf in the loop connecting the detector to the experiment. See Figure 1(b). This can be cured by using magnetic shielding, or by using twisted pair cables. An example of inductive coupling is the noise due to 50 Hz main power line, called “line interference” or the 50 Hz hum. It can be reduced by placing a notch filter which is centered around the line frequency.

Resistive Coupling Currents through common ground connections can give rise to noise voltages. This arises when we have two different grounding points which are not at exactly

the same potential. Thus, the detector measures the voltage across the experiment plus the voltage due to noise current passing through the finite resistance of the ground bus. This is shown by V_G in Figure 1(c). Resistive coupling can be cured by grounding everything to the same physical point.

Prudent laboratory practices can reduce noise sources to a manageable level and the lock-in amplifier employing the technique of phase sensitive detection can be used to recover signals which may still be buried in noise. We now discuss the working of a simple lock-in amplifier.

3.2 Internal working of a Lock-in Amplifier

The lock-in amplifier is used to detect and measure very small AC signals. A lock-in amplifier can make accurate measurements of small signals even when the signals are obscured by noise sources which may be thousand of times larger. Essentially, a lock-in is a filter with an extremely narrow bandwidth which is tuned to the frequency of the signal. Such a filter will reject most unwanted noise to allow the signal of specific frequency to be measured. For complete details the reader is referred to [3]. Figure 2 shows the approximate internal schematic of a typical lock-in amplifier.




Figure 2: Inner schematic of a lock-in amplifier. PSD represents the phase sensitive detector.

The lock-in amplifier is basically a phase sensitive detector (PSD) which comprises a mixer followed by a low pass filter. These components are explained here.

3.2.1 Mixer

The lock-in technique requires a reference frequency. We synchronize or derive the signal of interest from a suitable reference signal. The PSD operates by multiplying two signals together.

The demodulation process can be analyzed in two ways: graphically as well as mathematically. In the graphical method, illustrated in Figure (3) different possibilities for the signal are shown. In (a), (b) and (c), the signal and the reference are in phase, 90° out of phase and 0° out of phase respectively. First, let's consider case (a). the mixer output is positive. In the




Figure 3: The output V_o is determined by multiplying the signal and the reference wave. The subfigure (a), (b) and (c) show the signal and the reference in phase, 90° out of phase and 0° out of phase respectively.

negative half cycle of the reference, the mixer output is again positive. Therefore, when the reference and sinusoids are in phase, the mixer output is a full wave rectified sinusoid, whose dc component is proportional to the input signal.

Q 1. Explain the output when the input signal is 90° and 180° out of phase with the reference signal. What are the average values of the mixed signals?

We can also describe the mixing operation mathematically. Consider a sinusoidal input signal,

$$V_{in} = A \sin \omega t \quad \text{with} \quad \omega = 2\pi f, \quad (1)$$

where f and ω are the angular frequencies of the signal. For simplicity, we assume a reference signal of amplitude B , frequency Ω at phase ϕ relative to the input signal,

$$V_{ref} = B \sin(\Omega t + \phi). \quad (2)$$

The mixer operates by multiplying the two signals together so, the output V_o will be,

$$V_o = A \sin(\omega t) B \sin(\Omega t + \phi) \quad (3)$$

$$= \frac{AB}{2} (\cos((\omega - \Omega)t + \phi) - \cos((\omega + \Omega)t + \phi)), \quad (4)$$

showing that the mixer output comprises two AC signals, one at the difference frequency $(\omega - \Omega)$ and other at the sum frequency $(\omega + \Omega)$. If the reference frequency is equal to the frequency of input signal i.e. $\omega = \Omega$, a sinusoidal output is obtained with some DC offset Figure (3).

$$V_o = \frac{AB}{2} (\cos(\phi) - \cos(2\Omega t + \phi)). \quad (5)$$

So, the output V_o is proportional to the magnitude of input signal A , the cosine of angle between input and reference and is modulated at twice the reference frequency.

Q 2. Suppose the reference signal is a square wave whose amplitude is B and average value is $B/2$. What is the output from the mixer in this case?

3.2.2 Low Pass Filter

If the output voltage from the mixer, V_o is passed through a low pass filter whose cut off is $\ll 2\Omega$, the sinusoidal component is removed and we are left with the DC part only,

$$V_o = \frac{AB}{2} \cos(\phi). \quad (6)$$

Noise close to reference frequency contributes to a small magnitude to mixer output. The noise rejection depends upon the low pass filter bandwidth and the roll-off. As the bandwidth is made narrower, noise rejection is improved.

4 Finding an unknown resistance

In this section, you will determine the unknown tiny resistance R of a piece of jumper wire. This will be achieved by connecting it in series with a known resistance of some large value, say 10 K. With your understanding of the lock-in amplifier, find the resistance R . Compare your results with a nominal estimate obtained from a resistivity of $1.1 \times 10^{-7} \Omega\text{m}$. We expect you to draw a graph between two variables of your choice, and estimate R from the graph. We also expect a value for the uncertainty in R .

5 Finding a weak optical signal buried in noise (optional)

5.1 Objective

In this optional section, students will detect a weak optical signal using an optical chopper and a lock-in amplifier. When two polarizers are placed nearly crossed with each other, we obtain an extremely small signal on the photodetector which cannot be measured by an ordinary voltmeter. The lock-in amplifier can be used for the detection of such small signals, which are close to, but not quite zero. This extremely small signal will otherwise be washed out by ambient light.

The signal coming from the continuous wave laser is a DC signal whereas we require a modulated signal for the lock-in to operate. For this purpose, we use an optical chopper.

5.2 Optical Chopper

Modulation is the key step that enables the use of the lock-in amplifier. Modulation is the conversion of the DC signal to an AC signal with a defined frequency. In the present experiment, an optical chopper is used to square wave modulate the intensity of an optical signal. It comprises a chopping wheel, a motor and a speed control mechanism. The chopping wheel/blade is a rotating metallic disk with slots cut in the blade. The speed and number of slots determine the modulating frequency of the chopper. The modulation not only differentiates against noise but also discriminates against background light of constant intensity. Our optical chopper also provides an output frequency equal to the modulation frequency, which serves as the reference frequency for the lock-in amplifier.

5.3 Procedure

Carry out the following procedure to measure a weak signal (obtained by crossing two polarizers), using the lock-in amplifier and optical chopper.




Figure 4: Measurement of a weak optical signal using an optical chopper and lock-in amplifier.

1. Connect the setup as shown in Figure 4.
2. Switch on the optical chopper. Set the frequency to around 230 Hz. Keep the polarizers nearly crossed. The reference signal for the lock-in amplifier is provided by the chopper.

Precaution Always check the “OVL D” blinking on the lock-in amplifier. If it does, adjust the ‘sensitivity’ button.

Q 3. Connect the photodiode’s output to the provided I/V convertor. Observe the I/V’s output and the reference signal on the oscilloscope operating in the ‘dual mode’. Describe your observations.

Q 4. Adjust the polarizers to bring about a peak output voltage of 100 mV from the photodiode, as observed on the oscilloscope. Then, connect the signal to the lock-in amplifier.

Adjust the sensitivity, pre and post filters time constants and find the output. The lock-in amplifier shows an rms value in volts.

Q 5. Find the Fourier series of the I/V's output of Q2 and find the coefficient of the fundamental frequency. Correlate it with the measurement on the lock-in amplifier of Q3.

Q 6. Calculate the current being produced by the photodiode using the result of Q3.

Q 7. Now, connect the photodiode's output to the current input of the lock-in amplifier. Use the 10^6 V/A conversion factor to find the current. Verify the two readings.

Q 8. Now, cross the polarizers rendering the output on the oscilloscope immeasurable. Then, connect the output to the lock-in amplifier and see if you can still measure the output voltage. This step should be the most revealing aspect of the lock-in amplifier: its ability to detect a feeble voltage which is otherwise immeasurably small. Fully describe your observations and the settings of the lock-in amplifier.

Q 9. Find the output of the lock-in amplifier when the **REFERENCE MODE** is adjusted to $2f$. What does the output represent?

Q 10. In Q7, what is the optical power of the laser light after passing through the crossed polarizer?

Q 11. Find the output voltage using lock-in amplifier at the chopping frequency of 100 Hz. Can you explain the cause of the strange reading at this frequency?