

# Dynamics of water discharge

investigated with PhysLogger\*

Umar Hassan, Amrozia Shaheen, Mohsina Asif, and Muhammad Sabieh Anwar

LUMS Syed Baber Ali School of Science and Engineering

August 2, 2021  
Version 2021-v1

This experiment demonstrates the dynamics of water flowing out of a tank. Using the rate change of mass for an emptying cylinder, we investigate the application of Bernoulli's equation and the resulting Torricelli's theorem. We also observe the effects of constriction on the parcel of water flowing out of the tank. Students will investigate fluid dynamics, pressure and will relish how a phenomena as simple as water flowing out from a tank can lead to rich dynamics that can be explored mathematically.

## KEYWORDS

Pressure · viscosity · Bernoulli's equation · Torricelli's law · Laminar Flow · Continuity Equation

## 1 Conceptual Objectives

In this experiment, we will,

1. understand and apply Bernoulli's equation,
2. understand Torricelli's law,
3. understand the continuity equation,
4. learn how to numerically differentiate data, and
5. make plots of variables derived from directly measured quantities.

---

\*No part of this document can be re-used without explicit permission from Dr. Muhammad Sabieh Anwar.

## 2 Theoretical background

A fluid is a collection of molecules held together by weak cohesive forces. Usually liquids and gases are termed as fluids because they deform in response to external forces. Some general properties of fluid flow are summarized here.

1. **Steady or non-steady:** The flow of a fluid is described by pressure, density and flow velocity at every point of the fluid. If these variables are constant in time then the flow is steady.
2. **Compressible or incompressible:** If the density of a fluid remains constant and does not depend on  $x, y, z$  and  $t$ , then the flow is incompressible.
3. **Viscous or non-viscous:** Viscosity is the resistance towards flow. When a fluid flows such that there is no energy dissipation, then it is non-viscous flow. Such a flow is really an idealization.
4. **Rotational or irrotational:** If any element of the fluid does not rotate about an axis through the center of mass of the element, then the flow is irrotational.

### 2.1 Pressure inside a fluid

Consider a small segment of the fluid of density  $\rho$  at a distance  $y$  above some reference level as shown in Figure 1(a). This segment is a thin disk with thickness  $dy$  and area  $A$ , as illustrated in part (b) of the diagram. The mass of the element is  $dm = \rho dV = \rho A dy$  and its weight  $W = (dm)g = \rho g A dy$ . Since there is no acceleration the net vertical force is zero,




Figure 1: A static fluid. (a) Small element at rest inside the fluid, (b) forces acting on a small element.

$$\Sigma F_y = PA - (P + dP)A - \rho g A dy = 0, \quad (1)$$

yielding,

$$\frac{dP}{dy} = -\rho g. \quad (2)$$

This equation describes the variation of pressure with elevation above some reference level.


As the height increases ( $dy$  positive), the pressure decreases ( $dP$  negative). For an incompressible and homogeneous liquid with difference in height, the pressure difference is found by integrating the Equation (2)

$$P_2 - P_1 = -\rho g(y_2 - y_1), \quad (3)$$

and if the liquid has a free surface exposed to the atmospheric pressure  $P_o$ , then,

$$\begin{aligned} P_o - P &= -\rho g(y_2 - y_1), \\ P &= P_o + \rho gh, \end{aligned} \quad (4)$$

where,  $y_2 - y_1 = h$ . This shows that the pressure in a liquid increases with depth but would be same at all those points that are on the same level.



## 2.2 Bernoulli's equation

In our discussion on pressure, we have seen how pressure depends on the weight of the fluid above a level. However, pressure will also change with speed and elevation. You must have noticed how an object could be pulled into the wake of a fast-moving train; or by narrowing the hose of a water pipe, the stream of water can go further.

**Q 1.** Why do hordes of birds fly in a characteristic V-shaped pattern?

When a fluid moves through a region in which either the speed of the fluid or elevation above the earth's surface changes, the impact is that the pressure in the fluid changes. The relationship between fluid speed, pressure and elevation was first derived by Daniel Bernoulli in 1738. Bernoulli's equation, a fundamental relation in fluid mechanics is derivable from basic laws of Newtonian mechanics, as well as from the work-energy principle which stems from the conservation of energy.

Consider a steady, incompressible and nonviscous flow of a fluid through a pipeline from the position shown in Figure 2(a) to (b). The portion at the left has a cross sectional area  $A_1$  and at an elevation  $y_1$  from some reference level. A mass of fluid  $\Delta m$  gradually rises and after time  $\Delta t$ , it moves to the right end with cross sectional  $A_2$ , at an elevation  $y_2$ .

According to the work-energy theorem, the work done by the resultant force that acts on a system is equal to the change in kinetic energy. Assuming that there is no viscous force, the only forces that do work on the system are the pressure forces and the force of gravity. The net work done on the system by all the forces is,

$$W = P_1 A_1 \Delta l_1 - P_2 A_2 \Delta l_2 - (\Delta m) g(y_2 - y_1). \quad (5)$$




Figure 2: A fluid is flowing through a pipe from position (a) to (b). The net effect is the transfer of the element ( $\Delta m$ ) from the left to the right end. We calculate the work done in this transfer process.

This is the work done as the mass  $\Delta m$  displaces from (a) to (b). The pressure force  $P_2A_2\Delta l_2$  bears a negative sign because its direction is opposite to the horizontal displacement  $\Delta l_2$ . The gravitational force is also negative because it acts in a direction opposite to the vertical displacement. As  $A_1\Delta l_1 = A_2\Delta l_2$  is the volume of the fluid ( $\Delta V$ ) displaced, we can replace this  $\Delta m/\rho$ . The change in kinetic energy, therefore is,

$$\begin{aligned}
 \Delta K &= \frac{1}{2}\Delta m v_2^2 - \frac{1}{2}\Delta m v_1^2 \\
 &= \frac{1}{2}\Delta m (v_2^2 - v_1^2) \\
 &= P_1A_1\Delta l_1 - P_2A_2\Delta l_2 - (\Delta m)g(y_2 - y_1)
 \end{aligned} \tag{6}$$

This can be rearranged to give,

$$\frac{1}{2}\Delta m(v_2^2) + P_2 A_2 \Delta l_2 + (\Delta m)gy_2 = \frac{1}{2}\Delta m(v_1^2) + P_1 A_1 \Delta l_1 + (\Delta m)gy_1 \quad (7)$$

Dividing each side by the respective volume of the element  $A_2(\Delta l_2) = A_1(\Delta l_1)$ ,

$$P_2 + \frac{1}{2}\rho v_2^2 + \rho gy_2 = P_1 + \frac{1}{2}\rho v_1^2 + \rho gy_1. \quad (8)$$

The above equation is often expressed as,

$$P + \frac{1}{2}\rho v^2 + \rho gy = \text{constant}. \quad (9)$$

This is a statement of Bernoulli's equation.

The relation in Equation (8) can be modified in many different ways depending upon the situation. This leads to interesting corollaries. If the fluid is at rest i.e.  $v_2 = v_1 = 0$  then,

$$P_1 + \rho gh_1 = P_2 + \rho gh_2, \quad (10)$$

where the term  $(P + \rho gy)$  is called the *static pressure*.

Likewise, if both ends of the pipe are placed at same height then Equation (9) can be re-written as,

$$P_1 + \frac{1}{2}\rho v_1^2 = P_2 + \frac{1}{2}\rho v_2^2, \quad (11)$$

showing that high speeds corresponds to low pressures. The term  $\frac{1}{2}\rho v^2$  has dimensions of pressure and is called *dynamic pressure*.

## 2.3 Continuity equation

The equation of continuity for incompressible fluids states that,

$$\rho Av = \text{constant} \quad (12)$$

where  $v$  is the velocity and  $\rho$  is density of the fluid. This relation is easy to understand. Consider Figure (3) which shows a tapered horizontal pipe. The area at the left end is  $A_1$  and at the right is  $A_2 < A_1$ . In a unit time  $\Delta t$ , a mass of liquid  $\Delta m$  is transported between the ends. Since the fluid cannot be compressed, we must conserve the mass of fluid transferred, otherwise the liquid will turn denser in some regions and rarer in others. Therefore,

$$\rho A_1 \Delta l_1 = \rho A_2 \Delta l_2, \quad (13)$$

which dividing by  $\Delta t$  yields the equation of continuity, (12).




Figure 3: A tapered horizontal pipe. The horizontal velocity vectors are depicted by arrows.

**Q 2.** A giraffe needs a strong heart because of its long neck. Suppose the difference of height between the aortic valve (the place where the arterial blood comes out of the heart) and the head of a giraffe is 2.50 m, and the artery leading from near the aortic valve to the head has constant cross section all the way to the head. Blood is an incompressible fluid with density  $1.0 \text{ g/cm}^3$ . Assume the pressure at the head is zero.



- (a) What is the minimum required pressure at the aortic valve? Compare this pressure to the peak output pressure of the human heart ( $1.6 \times 10^4 \text{ Pa}$ )?
- (b) What would be the effect on the giraffe if the artery diameter narrowed down as it approached the brain?

## 2.4 Water discharge from a cylinder

A cylinder contains water that flows out from a narrow circular orifice at a fixed height  $y_2$  from the base. The orifice has a small area  $A_2$  compared to the cross sectional area  $A_1$  of the cylinder. As time progresses, the level of the water  $y_1(t)$  in the cylinder descends and water issues out with a speed  $v_2(t)$ . Let's apply Bernoulli's law to points **1** and **2**. Note that at the orifice, the jet of water is also exposed to the atmospheric pressure  $P_o$ . From Bernoulli's principle,

Consider  
Fig (4)

$$P_o + \frac{1}{2}\rho v_1^2 + \rho g y_1 = P_o + \frac{1}{2}\rho v_2^2 + \rho g y_2 \quad (14)$$

Since  $A_1 \gg A_2$ ,  $v_1 \approx 0$ , leading to

$$\begin{aligned} \frac{1}{2}\rho v_2^2 &= \rho g(y_1 - y_2) \\ v_2^2 &= 2gh, \end{aligned} \quad (15)$$




Figure 4: A cylinder with water flowing out from a narrow orifice at a fixed height  $y_2$  from the base.

that shows the relationship between the speed  $v_2(t)$  and the instantaneous head  $h(t)$  of water *above* the orifice.

## 2.5 Torricelli's Law

Torricelli's Law describes the relationship between the velocity of fluid leaving the cylinder  $v_2$  and the height  $h$  of the fluid. This relationship is given in Equation (15). In its simplest form, the speed,  $v$ , of a liquid flowing under the force of gravity out of an opening in a tank is proportional to the square root of the vertical distance,  $h$ . The speed of efflux is independent of the direction of flow. The theorem is named after Evangelista Torricelli, who formulated it in 1643. Notice that this speed is identical to the speed acquired by a mass falling under gravity through a height  $h$ .

In the experiment, you will observe if a linear relationship between  $v_2^2$  and  $h$  exists. Furthermore,  $v_2^2$  will in fact be observed to be smaller than  $2gh$ . The discrepancy will be accounted for by water's viscosity, and the effective narrowing of the orifice.

## 3 The Experiment

### 3.1 Preparation

You are provided a graduated cylinder with an orifice at a fixed height  $y_2$  from the base. Place it on the provided PhysLoad balance which is connected to the PhysLogger. Connect the PhysLogger to your PC's USB port and open the PhysLogger app. Using a USB-C

cable, connect the PhysLoad to one of the digital channels of the PhysLogger ( $P_A$ ,  $P_B$ ,  $P_C$ , or  $P_D$ ). Further, place the provided plastic box in line with the orifice to collect the water being discharged and running through the plastic water guide tube. Complete the assembly as shown in Figure 5.




Figure 5: The experimental assembly for observing the rate of discharge from a graduated cylinder.

In the PhysLogger app, select ‘PhysLogger’ and then click **Continue**. In **Options > Preferences**, define an appropriate destination folder where the experimental data will be stored. When prompted with ‘I want to’, select ‘Measure’. Next, scroll and select ‘Force’ and click on ‘Proceed’. A signal should appear on the plot. **Refer to Figure 6 and Table 1 for a description of the different relevant widgets and icons of the PhysLogger workspace.** Click on the ‘Time’ icon and choose a sampling rate of 0.2 Hz which means that a sample is acquired after every 5 s. Once these settings are completed, PhysLogger will begin plotting data acquired from PhysLoad.

Once you have set your axis you should tare the weight of graduated cylinder. For that purpose, click on the PhysLoad icon and select ‘Tare’. Alternatively, you may tare using the hardware **Tare** button on the PhysLogger. “Taring” means setting the reading of PhysLoad to zero. After you have your desired data acquired and displayed in PhysLogger’s window, you can export the data using the ‘Save’ option. (Quantities need to be dragged and dropped on the respective table’s columns before you click on save the table.)






Figure 6: The workspace of the PhysLogger Desktop App. Refer to Table 1 for a description of all icons and widgets.

No.	PhysLogger widgets and their descriptions
1	<b>Overlap layout</b> Displays overlapping graphs when more than one graph is being plotted.
2	<b>Cascade layout</b> Displays cascaded graphs when more than one graph is being plotted.
3	<b>Plot a new graph</b>
4	<b>Take screenshot</b>
5	<b>Edit layout</b>
6	<b>Time setup</b> Allows adjustment of sampling frequency.
7	<b>Analog input channels A–D</b> Indicate the voltage readings of the four analog input channels of PhysLogger.
8	<b>Output channels A and B</b> Indicate the voltage readings of the two output channels of the PhysLogger.
9	<b>PhysInstrument icon</b> Any PhysInstrument connected to the PhysLogger will appear here and allow the user to adjust instrument-specific settings (such as <b>taring</b> a PhysLoad).
10	<b>Clear live plot and recorded data</b>
11	<b>Pause data logging and live plotting</b>
12	<b>Record</b> To be used if data needs to be collected for only a specific, predefined period.
13	<b>Save</b> Saves and exports the plotted data points.
14	<b>Workspace setup</b> Allows the present hardware state of the PhysLogger to be saved for later use.
15	<b>Create a new page</b>
16	<b>Widgets and advanced tools</b>
17	<b>Plotted quantity label</b> Allows the units of the plotted quantity to be changed and the autoscaling of vertical axis to be toggled.
18	<b>Display options of the plot</b>
19	<b>Options</b> Allows the user to set <b>Preferences</b> and exit the PhysLogger application.

Table 1: Different widgets of PhysLogger (numbered as per Figure 6) and their functions.

## 3.2 Experimental procedure and analysis

In this experiment, you are required to verify a linear relationship between  $v_2^2(t)$  and  $h(t)$ . The balance will return the mass of water at each time step. Refer to Table 2 for the description of important values you will be needing in the analysis.

Apparatus Parameters	Values
Graduated cylinder diameter	$(69 \pm 1)$ mm
Diameter of the orifice	$(1.9 \pm 0.2)$ mm

Table 2: Various parameters pertaining the apparatus that will be required for analysis.

Write a computer programme that converts the obtained mass to the mass flow rate. From your programme, you should be able to compute the following quantities.

1. Rate of decrease of mass versus time.
2. Speed  $v_2$  of the discharging water jet versus time.
3. Height  $h$  of the descending water level versus time. Note that  $h$  is the head above the orifice and the orifice is at a height of 300 mL above the base of the cylinder.
4. Plot  $v_2$  versus  $h$ .
5. Plot  $v_2^2$  versus  $h$ . Fit this data to a linear relationship. Is it a good fit? Does the data corroborate Torricelli's theorem? Ideally the slope of the  $v_2^2$  versus  $h$  graph should have a slope  $2g$ . Is the slope of your data smaller or greater than  $2g$ ? How do you account for the difference?