



Phase Sensitive Faraday Rotation in TERBIUM GALLIUM GARNET crystal and various Diamagnetic liquid Samples

Supervisor: Dr. Saadat Anwar Siddiqi
Co-Supervisor: Dr. Muhammad Sabieh Anwar

Presented by: Aysha Aftab
Roll. No. 0702
M.Phil (2007-2009)



Outlines

- Magneto optics
- Polarization of light
- Jones Calculus
- Faraday Rotation
- Lock-in Amplifier ; PSD
- Why PSD in Faraday rotation?
- Schematic of the experiment.
- Results
- Determination of V using higher harmonics.
- References




Magneto Optics

Historical background

- 1813 Morrichini
- 1814 Faraday
- 1826 S. H. Christie
- 1834 Faraday
- 1844 “Magnetic force and light were proved to have relation to each other”.



Magneto Optics (cont...)




- Faraday rotation
- Kerr effect

Polarization of light


1. Linear polarization

- Orientation of **E** field remains constant, magnitude and sign varies.



Unpolarized light

- Plane of oscillation, containing **E** and **K**
- Electromagnetic wave with **E** field oscillating parallel to y axis.






Linearly polarized

Linear polarization (cont...)

Horizontally polarized light propagating in
z direction



$$E_x(z, t) = \hat{i} E_{ox} \cos(kz - \omega t)$$



Vertically polarized at phase difference ε ,


$$E_y(z, t) = \hat{j} E_{oy} \cos(kz - \omega t + \varepsilon)$$

Linear polarization (cont...)



Superposition of the two

$$\mathbf{E}(z, t) = \hat{i} E_{ox} \cos(kz - \omega t) + \hat{j} E_{oy} \cos(kz - \omega t + \varepsilon)$$





2. Circular polarization



- Two orthogonal waves have equal amplitudes.
- Relative phase shift of 90°
- Direction of \mathbf{E} is time varying, magnitude remains constant.

a) Right Circularly polarized light

- relative phase difference of $-90^\circ + 2m\pi$
- rotating clockwise

$$\mathbf{E} = E_o [\hat{i} \cos(kz - \omega t) + \hat{j} \sin(kz - \omega t)]$$

2. Circular polarization (cont..)





Right circular

- Left circularly polarized, phase shift of $90^\circ + 2m\pi$
- rotating anti clockwise

$$E = E_o [\hat{i} \cos(kz - \omega t) - \hat{j} \sin(kz - \omega t)]$$


2. Circular polarization (cont..)



Left circular

- Linearly polarized is sum of R.C.P and L.C.P

$$E(z,t) = 2E_o \hat{i} \cos(kz - \omega t)$$



Jones calculus


- R. Clark Jones in 1941
- For perfectly polarized light

$$E(z, t) = \hat{i} E_{ox} \cos(kz - \omega t) + \hat{j} E_{oy} \cos(kz - \omega t + \varepsilon)$$

$$E(z, t) = \begin{pmatrix} E_{ox} \\ E_{oy} e^{i\varepsilon} \end{pmatrix}$$

- For linearly polarized light

$$E(z, t) = \begin{pmatrix} \cos \alpha \\ \sin \alpha \end{pmatrix}$$





Jones calculus (cont...)


- Incident Jones vector E_i is related to transmitted Jones vector E_t through matrix, J

$$E_t = JE_i$$

- For beam passing through a series of optical elements

$$E_t = J_n \dots J_3 J_2 J_1 E_i$$

Faraday rotation



- Linearly polarized monochromatic light while transmitting through an optically inactive material, under the influence of an axial magnetic field, is rotated by an angle θ .
- Is non-reciprocal.



Faraday rotation (cont...)

- For uniform B field

$$\theta = VBd$$

- For non uniform B field

$$\theta = V \int_0^d B(z).dz$$

Where,

d= length of the sample


V= material parameter called Verdet constant,

- is a function of wave length of light.
- of the order of micro rad /G cm.


Faraday rotation (cont...)

Induced circular
birefringence

$$n_r \neq n_l$$



The plane of
polarization of
plane-polarized
light gets rotated





Phase Sensitive Detection (PSD)

$$V = 10\text{nV} \sin \text{ wave at } 20 \text{ kHz}$$

- Amplifier

$$\text{input noise} = 4 \text{ nV}/\sqrt{\text{Hz}}$$

$$Q = 1000 \quad \text{Bandwidth} = 1\text{MHz}$$

$$V_o = 1000 \times 10\text{nV} = 10\mu\text{V}$$

$$V_{\text{noise}} = 4\text{mV}$$



Phase Sensitive Detection (cont...)

- Band pass filter

at 50 kHz $Q = 1000$

$$\begin{aligned} V_{noise} &= \sqrt{500}(4\text{nV})(1000) \\ &= 89\mu\text{V} \end{aligned}$$

- PSD

Band width = 0.125 Hz


$$V_{noise} = 1.4\mu\text{V}$$



Lock in Amplifier

- A lock-in amplifier amplifies a small frequency band around a certain reference frequency.
- A lock-in can be used to
 - Measure sinusoidal voltage amplitudes and phase
 - Measure noise around a certain frequency.
- Consists of
 - Signal channel
 - Reference channel
 - PSD : Heart of lock-in amplifier



Lock in Amplifier (cont...)



Lock in amplifier, block diagram

Lock in Amplifier (cont...)

Mixer




$$V_{in} = A \sin(\omega t)$$

$$V_{ref} = B \sin(\Omega t + \varphi)$$

$$V_o = \frac{AB}{2} [\cos(\varphi) - \cos(2\omega t + \varphi)]$$

Lock in Amplifier (cont...)


Low pass filter




$$V_o = \frac{AB}{2} \cos(\varphi)$$

Noise close and closer to the reference frequency is removed.


Why PSD in Faraday rotation?



Why PSD in Faraday rotation? (cont...)




Noise and signal amplitude as a function of frequency.



Modulating the signal to a region of low noise.

Schematic of the experiment





- Light source
- Source of magnetic field
- Detecting devices

Schematic of the experiment (cont...)


Light source

■ Laser




(B) He-Ne gas laser

7/30/2010




(C) AlGaIn diode laser

Schematic of the experiment (cont...)



Linear polarizer





Malus's law

Source of magnetic field

- Dc source
large, expansive, bulky water cooled electromagnets.

Schematic of the experiment (cont...)





Schematic of the experiment (cont...)

- Ac source

Low magnetic field is required.

Signal of interest can be extracted successfully through PSD.

Solenoid, Helmholtz coil.

Helmholtz coil


Two identical coils with separation equal to their common radius.

$$B = B_o (1 + c_4 x^4 + c_6 x^6 + \dots)$$

In superposition

$$B = \left(\frac{4}{5}\right)^{3/2} \mu_o \frac{N}{a} i$$

Schematic of the experiment (cont...)



Magnetic lines of force of our Helmholtz coil drawn in Vizimag

Schematic of the experiment (cont...)




Schematic of the experiment (cont...)

■ Parameters of coil

- 18 gauge copper wire, $d=1.2$ mm
- $N=324$, $l=2.7$ cm , $R=1.5\Omega$
- $D1=6.5$ cm, $D2=10.2$ cm , $a=4.8$ cm
- $L=7$ mH

■ Connected a capacitor of $0.97\mu\text{F}$ to resonate the coil at 1.22 kHz to maximize the current and hence B

$$f = \frac{1}{2\pi} \sqrt{LC}$$



Schematic of the experiment (cont...)


Gauss meter





$$B = 6.7 \times 10^{-3} \text{ T}$$

$$Q = 38$$

Resonance in Helmholtz coil




Schematic of the experiment (cont...)



Magnetic field varies linearly with the current applied.

Schematic of the experiment (cont...)



7/30/2010

Setup for the experiment




Working principle

Jones vector of horizontally polarized light traveling in z direction

$$E_o = \begin{pmatrix} 1 \\ 0 \end{pmatrix} A_o \exp i(kz - \omega t)$$

For rotated polarized light, after sample

$$E_o = \begin{pmatrix} \cos \theta \\ \sin \theta \end{pmatrix} A_o \exp i(kz - \omega t)$$



Working principle (cont...)

After analyzer

$$E_o = \begin{pmatrix} \cos(\phi - \theta) \cos(\phi) \\ \cos(\phi - \theta) \sin(\phi) \end{pmatrix} A_o \exp i(kz - \omega t)$$

Intensity measured by photodiode




$$I = E^* . E$$

$$= A_o^2 \cos^2(\phi - \theta)$$


$$I \cong \frac{A_o^2}{2} (1 + 2\theta)$$

$$\cong I_0 + \Delta I \sin(\Omega t)$$

$$\theta = i_{ac} / 2i_{dc}$$




Resistance measurement through Lock-in Amplifier



$$R_a = \frac{V_c}{V_B} R_s$$




Using optical chopper with Lock-in Amplifier




Results



TGG at 405 nm





TGG at 405 nm





TGG at 633 nm





CS₂ at 405 nm





CS₂ at 633 nm





Methanol at 405 nm





Methanol at 633 nm





Water at 405 nm





Water at 633 nm





Isopropanol at 405 nm



Isopropanol at 633 nm



Determination of V using higher harmonics





Summary

sample	Wave length (nm)	Verdet constant (rad/G cm)
TGG	633	$(0.146 \pm 0.003) \times 10^{-3}$
TGG	405	$(1.70 \pm 0.02) \times 10^{-3}$
CS ₂	633	$(0.97 \pm 0.07) \times 10^{-5}$
CS ₂	405	$(0.50 \pm 0.01) \times 10^{-4}$
Methanol	633	$(0.24 \pm 0.02) \times 10^{-5}$
Methanol	405	$(1.20 \pm 0.05) \times 10^{-5}$



Summary (cont...)

sample	Wave length (nm)	Verdet constant (rad/G cm)
Water	633	$(0.30 \pm 0.02) \times 10^{-5}$
Water	405	$(1.50 \pm 0.08) \times 10^{-5}$
Isopropanol	633	$(3.55 \pm 0.02) \times 10^{-6}$
Isopropanol	405	$(2.2 \pm 0.06) \times 10^{-5}$



References

1. Aloke Jain, Jayant Kumar, Fumin Zhou and Lian Li, “ A simple experiment for determining Verdet constant using alternating magnetic fields” Am. J. Phys. 67, 714-717 (1999) .
2. Eugene Hetch and A. R. Ganesan, “Optics”, 4th edition, Pearson Education, Inc, India, 2008.
3. Frank L. Pedrotti and Peter Bandettini, “Faraday rotation in the undergraduate advanced laboratory”, Am. J. Phys. 58, 542-544 (1990).
4. Frank J. Loeffler, “A Faraday rotation experiment for the Undergraduate physics laboratory”, Am. J. Phys. 51, 661-663 (1983).



Reference

5. K. Turvey, "Determination of Verdet constant from combined ac and dc measurements, "Am. J. Phys. 64, 1561-1567 (1993).
6. <http://www.enzim.hu/~siza/cddemo/edemo0.htm>
7. Lock-in amplifier, user Manual, Stanford Research System, SR 510, <http://www.thinksrs.com>.

THANK YOU